

# EMERGENCY QUEST




---

Giovanni Clemente Monna  
Itzel Barreto  
Mlađo Milunović  
Stefano Codazzo

# SYSTEM CREATED FOR IMPROVING THE QUALITY OF LIVE OF :

People with  
**DEMENTIA / ALZHEIMER** & their  
**CAREGIVERS**



# MONITOR HEALTH.

In case of  
agitation



Play music  
Change Illumination

Relax the  
patient



# CARRYING OUT EVERYDAY ACTIVITIES

Movement in  
dangerous  
places



pre-recorded  
messages

Alert the  
patient to go  
back



THE OVEN




# ENHANCE SAFETY

Patient leaves  
a defined area



send  
notification  
to app

Notify the  
caregiver



**IMPROVE  
PATIENT LIFE**




by

smart  
house



monitoring  
bracelet



## INTERACTING

APP accessible  
for both users



- track patient movement
- store sensor data

## SENSING

GPS location



motion and  
agitation

## REASONING

Differentiate between



agitation  
problems



motion  
problems

## ACTING



pre-recorded  
messages



modify the  
environment

# ROOM STATION

---

## HARDWARE

Single board computer  
connected to the Internet  
via Wi-Fi

Motion sensor

Bracelet: heart rate and  
respiration sensor

Patient's smartphone to  
track his position

Active speaker

Controllable lights

## SOFTWARE



---

Station, sensors  
and actuators

process sensed data

process GPS patient's  
position

change the environment




# MAIN SERVER

---

## HARDWARE


computer  
connected to the  
Internet via LAN.

---



## SOFTWARE

store previous user data  
  
send notifications to the  
app  
  
receive and store  
appointments  
  
process notifications  
  
provide a map within the  
mobile app



## SOFTWARE

## MOBILE APP

## HARDWARE

### USER'S CELLPHONE

used to access to the system front-end from the mobile application.

create a new user to start using it



show patient's current position

set preferences

receive notification messages

patient history

set appointments with playing messages



**THANK YOU  
FOR YOUR ATTENTION.**

